



# Lung Disease Classification from Chest X-Ray Images

Using Deep Ensemble CNNs and a Python Streamlit Interface

**Eshan Puri**

M.E. – Artificial Intelligence | Project ID: 8025340013

# Project Overview

## ⚠ The Problem

Radiologists face an increasing workload, leading to fatigue and potential diagnostic errors. Manual analysis of Chest X-rays (CXRs) is time-consuming and subjective.

## 🎯 The Objective

To build an **end-to-end AI system** that provides:

**Automated "Second Opinion":** High-accuracy multi-class classification.

**Interactive Tool:** A user-friendly web interface for real-time diagnosis.

**Comprehensive Reporting:** Automated generation of medical-grade PDF reports.

## Target Classes

Bacterial Pneumonia    COVID-19    Tuberculosis  
Viral Pneumonia    Normal

## Key Technologies

Python    PyTorch    Streamlit    Pandas

# Methodology: Deep Ensemble Learning

---

We employed Transfer Learning using three diverse architectures pre-trained on ImageNet. The final predictor is a **Top-2 Soft-Voting Ensemble**.



## DenseNet121

Feature Reuse, 6.96M Params



## EfficientNet-B0

Compound Scaling, 4.01M Params




## ResNet50


Residual Learning, 23.5M Params

## Why Ensemble?

Combines the probabilities of the two best models (DenseNet + EfficientNet) to reduce variance and improve robustness.



# User Interface: Single Prediction



## Key Features

- **Real-time Inference:** Immediate feedback upon image upload.
- **Advanced Options:**
  - *Enhance Image (CLAHE):* For improving contrast in low-quality X-rays.
  - *Test-Time Augmentation (TTA):* Averages predictions of original and flipped images for reliability.
- **Visual Confidence:** Clear display of the predicted class (e.g., "Normal") and confidence score (99.20%).

**Observation:** The interface abstracts complex model parameters, presenting only clinically relevant data.

# User Interface: Batch Processing


## Simulating Clinical Workflow

Radiologists often need to screen multiple patients rapidly. The **Batch Prediction** module addresses this need.


**Bulk Upload:** Supports drag-and-drop for multiple files simultaneously.

**Tabular Results:** Displays Filename, Prediction, and Confidence in a structured table.

**Data Export:** Includes a "Download CSV" feature for offline record-keeping and analysis.



# User Interface: Evaluation & Metrics



## Transparency & Trust

Black-box AI is a barrier to medical adoption. This module provides transparency by exposing performance metrics directly in the app.

- **Dynamic Selection:** Users can select different models (DenseNet, EfficientNet, Ensemble) to compare performance.
- **Live Comparison Table:**  
Ensemble (Top-2): **0.855**  
EfficientNet-B0: **0.83**  
ResNet50: **0.81**
- **Visual Validation:** Displays Confusion Matrices and ROC curves dynamically.

# Error Analysis: Confusion Matrix

## Performance Breakdown


The confusion matrix for **EfficientNet-B0** reveals critical insights into the model's behavior.

### The "Viral-Bacterial" Challenge


The most significant error source is distinguishing between **Bacterial** and **Viral Pneumonia**.

- **45** Viral cases misclassified as Bacterial.
- **25** Bacterial cases misclassified as Viral.

*Insight: This mimics human expert difficulty, as these conditions share similar radiologic features (lobular consolidation vs interstitial patterns).*



## Discriminative Power: ROC Curve



## AUC Scores

The model exhibits near-perfect separability for distinct diseases.

Condition	AUC Score
Tuberculosis	1.00
COVID-19	1.00
Normal	0.99
Bacterial Pneumonia	0.96
Viral Pneumonia	0.94

# Automated Medical Reporting



## Standardized Documentation

The system auto-generates a PDF report for every prediction, closing the loop between diagnosis and documentation.

**Patient Metadata:** Auto-generated IDs and Timestamps.

**Quality Checks:** Includes *Blur Score*, *Brightness*, and *Contrast* to ensure the X-ray was readable.

**Disclaimer:** Explicitly states this is an AI-assisted analysis, not a final diagnosis.

AI-BASED DIAGNOSIS REPORT

ID: AI-1764180651

Prediction:

NORMAL

Confidence:

99.20%

Image Quality Metrics:

Blur Score:

116.21 (Clear)

Brightness:

115.50 (Optimal)

Contrast:

79.84 (Good)

Disclaimer: This is an AI-assisted analysis. Please consult a certified radiologist.

# Conclusion & Impact

## Summary

This project demonstrates the power of Python in creating a **holistic medical AI solution**. By moving beyond simple model

training to full-stack deployment, we bridge the gap between research and clinical utility.

## Future Scope

- **Explainability:** Integrate Grad-CAM heatmaps into the PDF report.
- **Feedback Loop:** Allow radiologists to correct predictions to retrain the model (Active Learning).
- **New Architectures:** Experiment with Vision Transformers (ViT).



85.4%

Final Ensemble Accuracy